

GE Healthcare

Probe Guide

Voluson™ Performance Series





4C-RS
H4000SR



RAB2-6-RS
H48681WR



IC9-RS
H48691PJ



RIC5-9A-RS
H48701EJ



12L-RS
H40402LY



3SC-RS
H45041DL

Description	Applications	Bandwidth	FOV	Availability
Curved – 2D				
Wideband Convex Probe	Abdomen, Obstetrics, Gynecology	2 – 5 MHz	81°	VP8, VP6
Curved – Real-time 4D				
Wideband Convex Ultra-light Volume Probe	Abdomen, Obstetrics, Gynecology	2 – 5 MHz	90° V 90° x 85°	VP8, VP6
Endocavity – 2D				
Wideband Microconvex Endocavitary Probe	Obstetrics, Gynecology, Transrectal	3 – 9 MHz	181°	VP8, VP6
Endocavity – Real-time 4D				
Wideband Microconvex Endocavitary Volume Probe	Obstetrics, Gynecology, Transrectal	4 – 9 MHz	180° V 180° x 120°	VP8, VP6
Linear – 2D				
Wideband Linear Probe	Small Parts, Peripheral Vascular, Pediatrics, Musculoskeletal, Breast	4 – 12 MHz	37 mm	VP8, VP6
Phased Array – 2D				
Wideband Phased Array Probe	Abdominal, Obstetrics, Cardiology, Pediatrics, Cephalic	1 – 4 MHz	90°	VP8

For probe care and cleaning information, visit www.gehealthcare.com/transducers.

Imagination at work

www.gehealthcare.com/transducers. Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information.

Data subject to change.

© 2016 General Electric Company.

GE, the GE Monogram and Voluson are trademarks of the General Electric Company.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

